| Model  | Go Go Elite Traveller® 2<br>Platinum (4-Wheel)  |
|--|---|
| Part Number                                      | SC442EMS  |
| Weight Capacity                                  | 300 lbs.  |
| Maximum Speed <sup>2</sup>                       | Up to 4.23 mph  |
| Ground Clearance <sup>4</sup>                    | 1.34" at motor<br>3.15" at mid-deck   |
| Turning Radius <sup>4</sup>                      | 37"   |
| Overall Length <sup>4</sup>                      | 39.37"  |
| Overall Width <sup>4</sup>                       | 19.45"  |
| Seat-to-Deck Height                              | 16.38 - 17.83"  |
| Seat Height                                      | 16.73"  |
| Tires  | Front: 2"x 7" Rear: 2.5"x 8"  |
| Range Per Charge <sup>2,3</sup>                  | 12 Ah: Up to 6.6 miles at 300 lbs.  |
| Product Weight                                   | 87 lbs. without batteries<br>106.5 lbs. with 12 Ah batteries<br>116.6 lbs. with 18 Ah batteries |
| Heaviest Piece<br>when Disassembled <sup>4</sup> | 39 lbs. (front section)   |
| Standard Seating                                 | Type: Mesh Back<br>Material: Black Vinyl<br>Dimensions: 17"x17"<br>Weight: 23 lbs.              |
| Motor Output                                     | Rear-wheel drive, sealed transaxle,<br>24-volt DC motor   |
| Brakes   | Regenerative and electromechanical  |
| Battery Requirements <sup>5</sup>                | (Two) 12-volt, sealed lead-acid,<br>non-spillable deep-cycle                                    |
| Battery Weight                                   | 12 Ah: 19.5 lbs. (standard)<br>18 Ah: 29.6 lbs.   |
| Battery Charger                                  | 2 amp off-board charger   |
| Warranty   | Lifetime limited on frame;<br>2-year limited on Electronics;<br>2-year limited on Drivetrain    |

# HCPCS Code: K0899

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.



(US) 800.800.8586 • (Canada) 888.570.1113

pridemobility.com







HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.





### FEATURE-RICH, TRAVEL-READY

The Go Go Elite Traveller® 2 Platinum is full of premium features to enhance your mobility experience. With 3-wheel scooter turning ability, magnetic charging, a convenient tiller mechanism, and feather-touch disassembly, the Platinum edition exudes luxury and performance mobility.

#### **Features**

- **EZ Turn technology** for the tight turning of a 3-wheel scooter and the stability of a 4-wheel (turning radius of 37")
- Charge 360 allows you to magnetically charge your scooter, requiring minimal dexterity and no need to line up any pins
- Angle-adjustable tiller mechanism can be operated with your hand or foot
- Feather-touch disassembly allows the unit to be easily separated into 5 lightweight pieces
- Front frame-mounted seat post for optimal stability
- Removeable footboard basket for additional storage
- 300-lb. weight capacity
- Maximum speed of up to 4.23 mph
- Front basket (standard)

2 sets of removable color-through shrouds included:

Cardinal Red (matte) Black (matte)

#### Accessories\*

- 2.1A XLR USB charger
- Cup holder
- Forearm crutch holder
- Military patches
- Oxygen holder
- Rear basket
- Rear mount 3" double cane/crutch holder
- Rear view mirror

- Saddle bag
- Safety flag
- Single rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- Additional shroud panels

## **Accessory Shroud Panel Colors**

(additional charge applies)



\* Only one rear accessory can be used at a time.